

Most Popular

	Verification*	Calibration	Calibration PLUS	Accredited
Step 1: Select One Service Level				
Manual Single Channel Pipette	\$30	\$43	\$49	\$74
Electronic Single, Repeater, Dispenser	\$31	\$49	\$54	\$77
Motorized Controller	\$31	\$49	N/A	N/A
6 Channel Pipette (manual/electronic)	\$51	\$51	\$57	\$172
8 Channel Pipette (manual/electronic)	\$51	\$74	\$80	\$221
12 Channel Pipette (manual/electronic)	\$51	\$86	\$92	\$284
16 Channel Pipette (manual/electronic)	\$51	\$97	\$103	\$347
24 Channel Pipette (manual/electronic)	\$51	\$109	\$115	\$438

Service Includes				
Leakage, Functional & Calibration Check	Included	Included	Included	Included
Preventive Maintenance (PM) **	N/A	Included	Included	Included
As Found Weighings	4x2	N/A	4x2	4x3
As Left Weighings	N/A	4x2	4x2	10x3
Calibration Report***	Included	Included	Included	ISO Format
Comprehensive CAL Label	Included	included	Included	Included

Step 2: Options				
3rd Point Weighings		<input type="checkbox"/> \$7	<input type="checkbox"/> \$14	Included
Extended 10 Weighings		<input type="checkbox"/> \$7	<input type="checkbox"/> \$14	<input type="checkbox"/> \$7
Multi-Channel Weighings (per channel)				
	<i>As Found</i>		<input type="checkbox"/> \$7	<input type="checkbox"/> \$14
	<i>As Left</i>	<input type="checkbox"/> \$7	<input type="checkbox"/> \$7	Included
• <i>3rd Point Weighings (per channel)</i>	<i>As Found</i>		<input type="checkbox"/> \$7	<input type="checkbox"/> \$14
	<i>As Left</i>	<input type="checkbox"/> \$7	<input type="checkbox"/> \$7	Included
• <i>Extended 10 Weighings (per channel)</i>	<i>As Found</i>		<input type="checkbox"/> \$7	<input type="checkbox"/> \$14
	<i>As Left</i>	<input type="checkbox"/> \$7	<input type="checkbox"/> \$7	Included
GMP/GLP Calibration Certificate		<input type="checkbox"/> \$5	<input type="checkbox"/> \$5	ISO Cert
DNA/RNA Decontamination (per pipette)		<input type="checkbox"/> \$15	<input type="checkbox"/> \$15	<input type="checkbox"/> \$15

* If pipette fails Verification, maintenance service will be performed under Calibration Plus level.
 ** PM (Preventive Maintenance) includes new Seal & O-ring as specified by manufacturer, piston polishing, thorough inside & outside cleaning.
 *** Included with returned pipettes / Uncontrolled copy of GMP/GLP available online at no charge. The laboratory implements the shared risk approach, based on the binary simple acceptance decision rule, as defined in ILAC G8 when making compliance statements.
 ****Calibration Warranty: 90 days labor / 1 year installed parts. Warranty not included with verification service.

Step 3: Calibration Cycle

MONTHS UNTIL NEXT CALIBRATION (Default is 6 Months)

3 Months 6 Months 1 12 Months Custom

PIPETTE QUANTITY

24 HOUR SERVICE (\$15 Per Pipette), limit of 25 pipettes per order, over 25 pipettes, please contact service@accuteklab.com. Transit time is not included, business days only.

I authorize repairs for up to \$100 per pipette. If repairs exceed \$100, please contact me for approval.

Paperless Account: Contact us for Qualer Access to Manage your pipettes and other assets for FREE!

Step 4: Decontamination Required

PIPETTES WERE DECONTAMINATED FROM RADIOACTIVE, CHEMICAL AND BIO HAZARDOUS MATERIALS WITH AN APPROPRIATE SOLUTION

CUSTOMER SIGNATURE

PLEASE PRINT THIS FORM AND INCLUDE EQUIPMENT LIST WITH YOUR SHIPMENT

Pipette Service Order Form

Step 5: Customer Information

FULL NAME

PHONE

E-MAIL

BILLING E-MAIL

SHIP TO ADDRESS

COMPANY NAME

DEPARTMENT

ADDRESS

CITY, STATE & ZIP CODE

BILL TO ADDRESS Check if same as "SHIP TO ADDRESS"

COMPANY NAME

DEPARTMENT

ADDRESS

CITY, STATE & ZIP CODE

Step 6: Payment Information

PO NUMBER OR CREDIT CARD NUMBER

EXP DATE /

Master Card Visa Amex

CVV Code Billing Zip Code

***Please include a PDF copy of your PO**

MY QUALER ACCOUNT (optional)

QUOTE WAS PROVIDED, (Please Complete Steps 3-7.)

Step 7: Shipping/Handling Options

USE MY CARRIER AND MY ACCOUNT # FOR SHIPPING

Carrier Account #

FREE 2-DAY FEDEX SHIPPING WITH QUANTITIES OF: 15 PLUS FOR CALIBRATION OR 20 PLUS FOR VERIFICATION for smaller quantities, shipping rates below apply:

Ground Shipping (one way) \$15

2nd Day Shipping (one way) \$25

Overnight Shipping (one way) \$45

Continental USA (other areas call for quote)

Rates shown above are based on Domestic Shipping Rates Only

Final Notes

- Include your equipment list: Print and attach it with your shipment.
- Complete and include this form with your shipment.
- For questions, contact us at service@accuteklab.com or 1-800-243-3232

Take Advantage of Our Trade-In Program: Receive Credit Toward Future Services or Products. Contact us at Service@accuteklab.com for the Trade-in Form and Instructions.

Flexible and Fast Calibration Services from Accutek Laboratories to Meet Your Lab's Needs.

- Customizable options are available, including multiple data points or weighings.
- 24-hour turnaround service for urgent needs (up to 25 pipettes per order).
- Batch shipments for larger orders with a 3-day turnaround time.
- Quick turnaround with two ISO 17025-accredited service labs across North America.

Our fast and efficient calibration options are designed to enhance your lab's productivity by minimizing downtime and ensuring precision. Please contact us at service@accuteklab.com or (800) 243-3232 to schedule your services.

Service Plans

All plans include:

- Full Preventive Maintenance (PM)
- Leakage, Functional & Calibration Check
- Adjustment to manufacturer's tolerances
- Calibration Certificate
- Calibration Label

As Found

As Returned

Preventive Maintenance

Calibration

Labs with basic service needs

Perfect for
Academic
R&D
Research Labs

N/A

4 Test weighings @
10, 100 %

Yes

Calibration Plus

When prior performance matters

Pharma
Biotech
Clinical
Research Labs

4 Test weighings @
10, 100 %

4 Test weighings @
10, 100 %

Yes

Accredited

For complete ISO 17025 documentation

GLP/GMP,
Regulated Labs

4 Test weighings @
10, 50, 100 %

10 Test weighings @
10, 50, 100 %

Yes

SHIP TO:



9477 Waples Street, Suite 120, San Diego, California 92121
www.Pipette.com | P: 1-800-243-3232 | F: 1-858-536-3342 | service@accuteklab.com

7690 First Place, Suite D, Oakwood Village, Ohio 44146

www.Pipette.com | P: 1-800-243-3232 | F: 1-858-536-3342 | service@accuteklab.com

Prices subject to change
QF11 020525 Pg. 2/2